

GENVEC

***“Comprehensive Characterization of the
293-ORF6 Cell Line***

Bryan T. Butman, GenVec, Inc.

GenVec, Inc.

- GenVec is a clinical stage biopharmaceutical company developing therapeutics for the treatment of major diseases, such as cancer, cardiovascular disease and macular degeneration.
- GenVec uses its patented adenovector technology to deliver medically relevant proteins, such as TNF-alpha, VEGF, and PEDF, directly to the site of disease.
- GenVec's adenovector technology is also being used to develop vaccines for infectious diseases.


GenVec's Pipeline

- Therapeutic product candidates:
 - Oncology: TNFerade for pancreatic and esophageal cancers - **TNF α**
 - Heart Disease: BIOBYPASS for severe coronary artery disease – **VEGF₁₂₁**
 - Ophthalmology: AdPEDF for wet age-related macular degeneration – **PEDF**

Vaccine Programs

- Collaboration with NIH/VRC for **HIV** and **SARS**
- Collaboration with U.S. Navy for **Malaria** and **Dengue Virus**


Multiple Product Opportunities




GenVec's Technology

- Replication-deficient Adenovirus vector technology
- “GV11” vector contains deletions in E1, E3 and E4 domains to assure replication incompetence
- Transgene inserted in place of E1 domain for protein production/delivery
- 293-ORF6 cell line developed to provide complementary E1 and E4 sequences

GV11(E1/E3/E4 Deleted) Vector



GenVec's 293-ORF6 Cell Line



- No overlapping sequences exist in the E4 regions of GV11 (E1/E4 deleted) vectors and 293-ORF6 cell line

293-ORF6 MCB Characterization

- Viral and microbial safety panel
- Cytogenetic analysis
- E1 and E4 copy number
- E1 and E4 chromosome localization
- PrP gene sequence
- PrP gene product protease sensitivity
- Tumorigenicity and Oncogenicity
- PCR identity assay

Safety Assessment of MCB

- **Microbial**: Mycoplasma, Sterility (B&F)
- **In Vitro**: Adventitious viruses (3 cell lines), Bovine Viruses (9 CFR), Porcine Parvovirus
- **In Vivo**: Inapparent viruses in animals
- **PCR**: AAV, CMV, EBV, HBV, HCV (RT-PCR), HHV-6/7/8, HIV-1/2, HPV, HTLV-I/II, Parvo B19, SV40, F-PERT (Retroviruses)
- **Biochemical**: Isoenzyme analysis (human origin)
- **Other**: TEM

Cell Aging Study

- “Aged Cell Bank” established to assess cell line stability
 - “P43” designation
 - 16 passages beyond the MCB (P27)
 - Represents theoretical expansion of WCB to 10,000L bioreactor stage
 - P43 bank cryopreserved for characterization studies

Characterization of MCB

- P27 (MCB) compared to P43 (aged cell bank) for the following characterization studies:
 - Cytogenetic analysis (ploidy)
 - Copy number of E1 and E4 insertions
 - Chromosome localization of E1 and E4 insertions
 - DNA sequence of PrP gene
 - Protease sensitivity of PrP gene product
 - PCR identity assay for 293-ORF6 cells

Cytogenetic Analysis

- P27 and P43 cells were expanded, treated with Colcemid, hypotonically lysed, fixed and then stained with Wright's reagent
- 100 cells in metaphase were examined for chromosome count
- Structural abnormalities and chromosome number were recorded
- Abnormalities classified as:
 - Chromatid-type aberrations
 - Chromosome-type aberrations
 - Severely damaged cells

Cytogenetic Analysis

Results:

	<u>P27</u>	<u>P43</u>
Median#/Metaphase	72	72
Range	66-80	54-78
Aberrations/cell	0.02	0.00

E1 and E4 Copy Number

- Q-PCR assay developed
 - E1 assay targeted Ad5 bases 3343 to 3411
 - E4 assay targeted Ad5 bases 34007 to 34073 (specific for ORF6 region)
 - Human GAPDH used as an index for cell number

E1 and E4 Copy Number

Results:

	<u>P27</u>	<u>P43</u>
E1 Copy #	6.5 ± 0.8	5.1 ± 0.7
E4 Copy #	2.6 ± 0.2	2.4 ± 0.4

E1 and E4 Localization

- FISH assay developed to detect E1 and E4 chromosome localization
 - Plasmids encoding full E1 region or the E4-ORF6 expression cassette were labeled with Texas Red 5'UTP by nick translation
 - P27 and P43 cells were expanded, treated with Colcemid, lysed, fixed and then stained with labeled probes and counterstained with DAPI
- Cells in metaphase were analyzed for the presence of red fluorescent signal
- Chromosomes were identified by DAPI

E1 and E4 Localization

Results:

P27

P43

E1 Local. 1 copy on #6
at 6q16-21

1 copy on #6
at 6q16-21

E4 Local. 2 copies on #19
at 19q13.3

2 copies on #19
at 19q13.3

PrP Gene Sequence

- Prion protein (PrP) open reading frame found on chromosome 20 was sequenced directly from total DNA
- PCR Primers:
 - AP-1 – Palmer, et al, 1996
 - AP-2 – Windl, et al, 1996
- Resultant amplicon was purified and sequenced

PrP Gene Sequence

- Expected 894 bp amplicon coding for 245 aa Prion protein, based on GenBank accession #U29185
- PrP-specific PCR products were amplified (918 bp) from both P27 and P43 cells
 - Sequences matched GenBank #U29185
 - Insertion of 24 consecutive bases between bp 339 and 340 in both sequences
 - Insertion represents known octapeptide repeat
- **There was no evidence of infectious PrP^{Sc}**

PrP Gene Product Protease Sensitivity

- Western blot assay employed that detects protease-resistant PrP gene product (Caughey, et.al, 1996, 1998)
- Assay performed on both P27 and P43 cells
 - Cells were cultured, PrP protein extracted, and subjected to Proteinase K digestion
 - Materials analyzed by Western blot using Mab 3F4 that recognizes both human and hamster PrP-res
 - Positive control = Hamster scrapies agent
 - Assay sensitivity $\leq 99\%$ PrP sensitive to protease

PrP Gene Product Protease Sensitivity

Results:

Controls – All spiked samples and controls showed presence of protease-resistant PrP

P27 – No evidence of protease-resistant PrP

P43 – No evidence of protease-resistant PrP

Tumorigenicity Study

- Goal: Determine whether addition of ORF6 plasmid increases tumorigenic potential of 293-ORF6 cells over parental 293 cells
- Animals: Athymic Nude (nu/nu) mice, 20 mice per group
- Test Article: 293 parental and 293-ORF6 (P43) cells injected at 10^7 , 10^5 , and 10^3 cells/animal
- Positive control: HeLa at 10^7 , 10^5 and 10^3 cells/animal
- Negative control: Syrian Hamster Embryo (SHE) at 10^7 cells/animal
- Study Duration: 140 days
- Endpoint: Histopathically-confirmed neoplasms

Tumorigenicity Study

Results:

HeLa Positive Control – Tumorigenic at 10^5 and 10^7 cells

SHE Negative Control – Not tumorigenic at 10^7 cells

293 Parent – Tumorigenic at 10^7 cells

293-ORF6 (P43) – Tumorigenic at 10^7 cells


Conclusion: No difference observed

Oncogenicity Study


- Purpose: To further assure safety of vaccines planned with 293-ORF6 cells
- Design:
 - Test article: P43 cells
 - Cell lysate
 - Purified DNA
 - Species: Newborn rats, newborn hamsters
 - Administration: Subcutaneous, equivalent of 10^7 cells/animal
 - Controls:
 - Negative – MRC-5
 - Positive – None
- Study duration: 5 months
- Endpoint: Histopathically-confirmed neoplasms

293ORF6 Identity Assay Design

	<u>E1</u>	<u>ORF6 EC</u>
293 cells	+	-
293/ORF6 (P27)	+	+
293-ORF6 (P43)	+	+
A549 cells	-	-
J1-20 cells	-	+



Serum Free Suspension 293-ORF6 Cell Stability With Ad_{Gv}EGR.TNF.11D



Summary & Conclusions

- Full characterization of GenVec's 293-ORF6 cell platform
 - Safety
 - Molecular characterization
 - Stability
- Submission of BMF to FDA
- Approval for use of 293-ORF6 cell line in the development of clinical vaccines

Acknowledgements

- Co-authors – **Alena Lizonova, Jason Gall, Perry Newton, Mike Sowers and Doug Brough**, GenVec; **Rebecca Sheets and Phil Gomez**, VRC/NIAID
- PCR analyses – **Audrey Chang**, BioReliance
- Cytogenetic analysis – **Carlos Sanchez**, Applied Genetics Laboratories
- E1 and E4 copy number and localization – **Elena Jouravleva**, Lark Technologies
- PrP sequencing – **Heath Knight**, Lark Technologies
- Protease resistance WB – **Byron Caughey**, Rocky Mountain Laboratories, NIAID
- Tumorigenicity – **Mary Whiteman and Dan Ragland**, BioReliance
- Cell identity – **Gerald Lohnas**, GenVec
- Project Management – **Victoria Haque**, GenVec

GenVec Products in Development

