

Fluorescent Labeling of Total RNA

Updated 01/28/04

1. Add primer to the RNA samples to be labeled.

- Make up fresh or thaw previously made Master Mixes B. (See Master Mix Table)

Prime the two RNA samples to be compared

RNA sample (30 µg total RNA per rxn) in dH ₂ O	15 µl
Oligo dT primer (1mg/ml)	1 µl
Total	16 µl

(See notes 1,2,6)

2. Heat and snap cool.

- Heat samples to 70°C for 5 min., then snap cool in an ice-water bath. Centrifuge briefly to collect the contents to the bottom of the tube.

3. Initiate two labeling reactions, one for Cy3 and one for Cy5.

Add these components in the order listed

RNase inhibitor diluted 2X	0.5 µl
Master Mix B	10.0 µl
dUTP-Cy3 OR dUTP Cy5	2.0 µl
Superscript II	1.5 µl

Total Volume should now be 30ul

(See Note 3)

- Mix the components well with the RNA/primer sample and centrifuge briefly to collect the contents to the bottom of the tube

4. Incubate reverse transcription reaction mixture at 42°C for 90 min

5. Degrade RNA

- After the reverse transcription reaction is complete, add 10µl of 1M NaOH and heat at 70°C for 15 min.
- Neutralize with 10µl 1M HCl.

6. Clean up labeled cDNA

- Combine the appropriate Cy3- and Cy5-labeled samples that will be co-hybridized onto the array.
- Add the combined probes to a Vivaspin 500 (10K cutoff) with 400 µl TE. Centrifuge at 12,000 X g (~10,000 rpm) for 8 min. Discard flow through.
- Repeat the wash step in the same Vivaspin filter at least three times.
- We recommend 2-3 additional washes in a fresh Vivaspin 500. Add 400 µl TE, centrifuge at 12,000 X g (~10,000 rpm) for 8 min. Discard flow through.
- Remove the sample retained in the filter chamber with a pipette and transfer to a new tube. Adjust the sample volume to 20 µl for the mouse and human standard microarray formats.

Notes

Note 1. This protocol is designed to produce two cDNA samples in parallel that are labeled with different fluorescent dyes – one with Cy3 and the other with Cy5. Total RNA is the template for the reverse transcription reaction that creates the cDNA copy and directly incorporates a fluorescent nucleotide (either Cy3-dUTP or Cy5-dUTP). These two differentially labeled cDNA pools are combined in Step 6 and co-purified to remove excess nucleotides, primers, salts and buffers. Adequate purification of the labeled cDNA is essential to reduce hybridization background fluorescence. We have found that labeling of 30µg RNA is appropriate for both Cy3 and Cy5.

Note 2. Keep the RNA and other reaction components on ice until instructed otherwise. Use RNase-free reagents and technique throughout the labeling procedure. The best results are typically generated by high-quality, pure RNA preparations. Users are encouraged to check the integrity and quantity of their RNA before labeling.

Note 3. Use of the RNase inhibitor is optional.

Master Mix B Table

Reagent	Stock concentration	Volume for 1X (Assumes 30µl reaction volume)	Bulk mixture for 100 reactions
5X RT buffer (Invitrogen)	5X	6 µl	600 µl
100 mM DTT	100 mM	3 µl	300 µl
dATP	100 mM	0.15 µl	15 µl
dCTP	100 mM	0.15 µl	15 µl
dGTP	100 mM	0.15 µl	15 µl
dTTP	100 mM	0.03 µl	3 µl
Water (RNase-free)		0.52 µl	52 µl

Master Mix B should be made up in sufficient volume to allow for accurate pipetting of each component. Master Mix B may be stored frozen. Distribute Master Mix B into aliquots to minimize the freeze-thaw cycles. Avoid thawing a Master Mix B tube more than twice after it is made. Thaw or create enough Master Mix solution in sufficient quantity to complete all reactions plus one

Note 4. Researchers must keep accurate notes on which samples are labeled with which fluorescent dyes. The table below defines some of the properties of each dye.

Cyanine dye	Cy3	Cy5
Visual color of dye	Pink	Blue
Excitation maximum	550 nm	650 nm
Emission maximum	568 nm	668 nm
Extinction coefficient ($M^{-1}cm^{-1}$)	150000	250000
Axon scanner channel	532 nm	635 nm
Standard pseudo colored display (by the software)	Green	Red
Vulnerabilities	Higher background than Cy5	Quenched by dense packing; more labile on dry slides

Appendix

Table A-1. Sources of materials and reagents for microarray labeling and hybridization

Vendor	Item	Catalog
Ambion www.ambion.com 800-888-8804	Superase In (RNase inhibitor)	2694
Amersham Pharmacia www.apbiotech.com 800-526-3593	Cy3-dUTP, 25nmol, quote 752BQ	PA53022
	Cy5-dUTP, 25nmol, quote 752BQ	PA55022
	dNTP set (dATP, dCTP, dGTP, dTTP) 4X25umol	27-2035-01
	Random hexamer pdN(6) (for prokaryotic RNA)	27-2166-01
Invitrogen , www.invitrogen.com , 888-584-8929, BPA# 263-00040377 ***	Superscript II reverse transcriptase, RNaseH-neg	18064014
	Oligo dT primer (12-20mer) 1 mg/ml	Poly.gf
	100mM dNTP set, 4X 25 μ mol	10297-018
Vivascience www.vivascience.com	Vivaspin 0.5ml Concentrator	VS0102

* Must reference quote # 752BQ for NIH discount from Amersham Pharmacia through 12/31/04.

** May require lead time, minimum order 5oz

*** BPA pricing is less than catalog listing in table

This table lists materials and products that will perform adequately in this procedure. However, this table does not represent an endorsement of any company or product. There are other products that can be substituted for the items listed with no adverse consequence on the results.