



Detection of Adventitious Viruses in Biologicals – A Rare Occurrence

Vaccine Cell Substrates 2004

June 30th, 2004

Raymond Nims

Rockville, Maryland



Detection of Viruses in Biologicals Using Screening Assays - Conceptual Basis

Principle: Viruses will replicate in a suitable host cell, expressing their presence through cytopathic effects, and/or hemagglutination or hemadsorption of specific types of erythrocytes

Indicator cell lines are selected on the basis of theoretical susceptibility to viruses of concern. In *in vitro* screening, a **human diploid** cell and a **primate** cell are used to detect viruses infectious for human cells, while a third indicator of a **similar species** as the cell substrate used in manufacture of the biological is used to detect viruses infectious for those cells.

In bovine and porcine virus testing, bovine and porcine indicator cells are employed.

- In USA: 9 CFR 113.53
- In EU: ICH Viral Safety Documents Q5D (1993), Q5A(1997)

- **In USA:** Driven by the Points to Consider Documents (guidelines formulated by the FDA, not laws)
 - Characterization of Cell Lines Used to Produce Biologicals (1993)
 - Manufacture & Testing of Monoclonal Antibody Products (1997)
 - Guidance for Human Somatic Cell Therapy & Gene Therapy (1998)
- **In EU:** ICH Viral Safety Documents Q5D (1993), Q5A(1997), European Pharmacopoeia Commission 1999 Tests for Extraneous Agents in Viral Vaccines for Human Use.

- Raw materials (serum, trypsin, other animal-derived materials) 9 CFR-compliant bovine and porcine virus screens
- Cell banks (WCB, MCB, EPC) In vitro virus screen, bovine or porcine virus screens
- Unprocessed bulk harvest (Lot Release) In vitro virus screen

- **Cellular vaccine** products may be evaluated in a manner similar to Cell Line Characterization
- **Live viral vaccines** (vaccinia, influenza, etc.) typically require neutralization with type-specific antisera prior to testing to avoid false positive effects in the screening assays.

Raw Materials:

- Serum (15% in medium)
- Trypsin (5 g, centrifuged, reconstituted in trypsin inhibitor)
- Other animal-derived materials (3x nominal concentration)

Cell banks: 10^7 cells/ml in conditioned medium

Bulk harvests: sample size (3 ml per indicator cell line) has no relationship to total harvest volume.

1. Host cells are exposed to test agent in flasks or plates
2. Cells are incubated for 2-4 weeks and observed under inverted microscope
3. Test cultures are compared to negative and positive control cultures
4. At termination, additional endpoints may be performed (hemadsorption/hemagglutination of erythrocytes, Immunofluorescence)

- Consideration of host range, cytopathic effect, hemadsorption/hemagglutination pattern
- Isolate remains infectious after **passage** to fresh cells?
- Conduct **PCR** testing for likely candidates
- Conduct **IFA** screening of infected cells
- Conduct **Electron microscopy** of infected cells

<u>Raw Material</u>	<u>Viruses Detected</u>
---------------------	-------------------------

Bovine serum	BVDV bovine polyoma virus
--------------	------------------------------

Porcine trypsin	PPV
-----------------	-----


<u>Biological Type</u>	<u>Viruses Detected</u>
------------------------	-------------------------

Cell Lines (MCB, WCB, EPC)	none
Gene Therapy Vectors	replication-competent adenovirus
Monoclonal Antibodies	none
Recombinant Proteins	
• non-CHO cell process	none
• CHO cell process	MMV, REO, Cache Valley Virus
Vaccines	none


Morphological Changes in 324K Cells Infected with REO Virus



REO-infected 324K cells
Day 20 10x magnification



Control 324K cells
Day 20 10x magnification



- packing of viral particles in crystalline arrays near nucleus at late stage of infection
- 65-75 nm diameter spherical particles -with double concentric ring structure
- 37000x magnification

- Raw materials testing can reduce potential for introduction of viruses in serum, trypsin, other animal-derived materials
- Adventitious virus screening is intended to detect gross contamination of raw materials, cell lines, and bulk harvests. Contaminants are rarely detected in cell substrates and bulk harvests. Information obtained from such occurrences is used to design cleaning and process (viral clearance) validation studies.