



Scope & Practicality of *In Vivo* Testing for Adventitious Agents

Vaccine Cell Substrates 2004

June 30, 2004


Mary D. Whiteman

Rockville, Maryland



- ◆ Oncogenes
- ◆ Potential To Stimulate Oncogenes
- ◆ Bacterial, Fungal and Mycoplasmal Contamination
- ◆ Viruses - Human
- ◆ Viruses - Rodent

Portals of Entry for Contaminants



Key Regulatory Guidelines and Documents

- ◆ United States
 - 21 CFR Parts 58, 200-299, 600-680
 - “Points to Consider” - CBER (1997)
- ◆ European Union
 - “Notes for Guidance” - CPMP (1998)
- ◆ International Conference on Harmonization (ICH)
 - ICH Viral Safety Document (1997)

Product**Contaminant**

Yellow Fever Vaccine

Avian Leukosis Virus, Hepatitis
B Virus

Poliovirus Vaccine

SV40

Poliovirus Vaccine

Live Poliovirus

Adenovirus Vaccine

SV40

Blood Products

HIV

Human Growth Hormone

Creutzfeldt-Jakob Disease

- Biopharmaceutical Products
 - Monoclonal Antibodies – Mouse, Rat, Human
 - Proteins from Genetically Engineered Mammalian Cell Lines (e.g. CHO)
- Gene Therapy Products (Viral Vectors)
- Human Blood Derived Products
- Biological Pharmaceuticals
 - Heparin (Bovine Origin)
 - Collagen

- Vaccines
- Xenotransplants
- Supplements/Raw Materials used in Pharmaceuticals
 - Bovine Serum
 - Sheep Blood
 - Peptones/Amino Acids
 - Bovine Insulin/Transferrin

1. Selection and testing of the source materials for the absence of virus.
2. Testing the product at appropriate stages of production for the absence of virus.
3. Testing the capacity of the production processes to remove or inactivate viruses.

- ◆ No single test can demonstrate the presence of all known viruses.
- ◆ All test systems require a minimum level of viral contamination to record a positive.


- If the sample contains a low concentration of virus and only a fraction of the sample is assayed, the fraction may test negative due to random and unequal distribution throughout the sample.

- Viruses that do not cause cytopathic or other noticeable effects in a cell culture system may be detectable in an animal system.

- 21 CFR Part 630.35 (Revoked)
- 21 CFR Section 630.16 (Revoked)
- Jacobs JP, McGrath DI, Garrett AJ and Schild GC. 1981. Guidelines for the acceptability, management and testing of serially propagated human diploid cells for the production of live virus vaccine for use in man. *J. Biol. Stand.*, **9**, 331-342
- European Pharmacopoeia Commission 1999. Tests for Extraneous Agents in Viral Vaccines for Human Use (General method N° 2.6.16.). European Pharmacopoeia 3rd edition Supplement: 50-51.

- Animal systems
 - Suckling mice
 - Adult mice
 - Guinea pigs
 - Embryonated hen's eggs

- **Experimental Design – Suckling Mice**



- Each pup injected
 - intraperitoneally (i.p.)
 - intracranially (i.c.)
 - per os (p.o.)
- Animals observed for 14 days.
- On day 14, a single pool of emulsified tissue of all surviving mice is prepared and passaged.
- Passage animals observed for 14 days.

- **Agents for which suckling mice are an efficient isolation system:**
 - Arboviruses
 - Coxsackie A viruses
 - Coxsackie B viruses
 - Herpes Simplex (type 1 and 2)
 - Rhabdoviruses (including rabies)
 - Togaviruses (LDV)
 - Junin
 - Herpes B
- **Agents for which suckling mice are of secondary efficiency as an isolation system:**
 - LCM
 - Lassa
 - Hantaan
 - Ebola
 - Vaccinia

- **Experimental Design – Adult Mice**

Day 0
Injection




Day 28
Sacrifice

- Each mouse injected
 - intraperitoneally (i.p.)
 - intracranially (i.c.)
 - per os (p.o.)
 - intranasally (i.n.)
- Animals observed for 28 days.

- **Agents for which adult mice are an efficient isolation system:**
 - Rhabdoviruses (including rabies)
 - Togaviruses (LDV)
 - LCM
- **Agents for which adult mice are of secondary efficiency as an isolation system:**
 - Arboviruses
 - Herpes Simplex (type 1 and 2)
 - Lassa


- **Experimental Design – Guinea pigs**




- Each guinea pig injected
 - intraperitoneally (i.p.)
 - intracranially (i.c.)
- Injection sites observed weekly for lesions.
- Animals observed for at least 28 days.
- For vaccine tests, animals observed for 42 days.

- **Agents for which guinea pigs are an efficient isolation system:**
 - Rhabdoviruses (including rabies)
 - LCM
 - Lassa
 - Junin
 - Marburg
 - Ebola
- **Agents for which guinea pigs are of secondary efficiency as an isolation system:**
 - Arboviruses
 - Vaccinia

- **Experimental Design – Allantoic Route**



- Eggs injected and incubated for 3 days.
- Allantoic fluids are tested for hemagglutinins.
- Fluids are pooled and passaged.
- Fluids from passage eggs are tested for hemagglutinins.



- **Experimental Design – Yolk Sac Route**


Day 0
Injection



Day 9
Viability
Passage



Day 9
Viability



- Each egg injected and incubated for at least 9 days.
- Examined for viability.
- Yolk sac material is pooled and passaged.
- Eggs incubated for at least 9 days.
- Examined for viability.

Isolation and Detection of Viruses – Embryonated Hen's Eggs

- **Agents for which hen's eggs are an efficient isolation system:**
 - Herpes Simplex (type 1 and 2)
 - Rhabdoviruses (including rabies)
 - Herpes B
 - Mumps
 - Influenza
 - Parainfluenza (types 1, 2 and 3)
 - Vaccinia
- **Agents for which hen's eggs are of secondary efficiency as an isolation system:**
 - Arboviruses

- Evaluation of Test Results
 - The test material will be considered not contaminated with adventitious viral contaminants if the following are met:
 - 80% of animals/eggs
 - Remain healthy
 - Survive the observation period
 - No evidence of viral infection

- Suckling mice
 - Death by natural causes
 - Runting
 - Inadequate maternal care
 - Cannibalization
 - Injection trauma
- Adult Mice/Guinea pigs
 - Death by natural causes/injection trauma
 - Neurological clinical signs due to i.c. injection
- Eggs
 - Hardiness (death with no assignable cause)
 - Bacterial contamination

- Gene Therapy Vectors
 - Toxicity
 - Infection
 - Hemagglutination
- Bulk Harvest Material
 - Media/Component toxicity
- Vaccine Viral Stocks and/or substrates
 - Substrate bioburden
 - Virus infects test system
 - Virus hemagglutinates
 - Viral neurovirulence
 - Antiserum toxicity

- Perform feasibility and/or qualification studies
 - Dilution
 - Neutralization
 - Antibiotics
 - Antiserum concentration
 - Antiserum control animals
 - Elimination of one or more test systems or routes of injection
 - Determination of viral breakthrough

Data from assays performed at Rockville, MD from
March 28, 2003 through August 31, 2003

Test System	Number of Assays Examined	Number of System Suitability Problems	System Suitability
Suckling Mice	98	7	93%
Adult Mice	98	1	99%
Guinea pigs	88	2	98%
Allantoic route	100	3	97%
Yolk sac route	100	5	95%

- Tests for Adventitious Contaminants
- Antibody Production Tests
- Tumorigenicity Test
- Rabbit Pyrogen Tests (USP and EP compliant)
- General Safety Tests (USP and EP compliant)
- TSE Clearance Bioassay

- Sarah Gilliam, Ph.D., BioReliance, Scotland
- Ramakrishna Velicheti, Ph.D., BioReliance, Rockville, MD
- Cynthia I. Sigler, B.A., BioReliance, Rockville, MD